ProDrive Standard/Mini

Highspeed Handpieces
Turbines dentaires
Dentalturbinen

PRODRIVE
SYSTEMS

Handpiece Performance Through Innovation

English / Français / Deutsch

Operating Instructions
Notice d'utilisation
Gebrauchsanweisung

ProDrive Standard/Mini Highspeed Handpieces – Operating Instructions

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Highspeed Handpieces

PRODRIVE
SYSTEMS

Handpiece Performance Through Innovation

English

Operating Instructions

We reserve the right to make any alterations due to technical improvements.

Caution: Federal law restricts this device to sale by or on the order of a Dentist

Dear customer,

ProDrive highspeed handpieces are designed for dental treatment performed at dental treatment centers and for laboratory use.

Please follow the warning and safety instructions provided to prevent personal injury and/or material damage. These instructions are individually marked WARNING, CAUTION or NOTE.

List of Contents

1	Technical Description	6	2	Operation	11
1.1	ProDrive Standard and Mini	6	2.1	Handpiece with ProDrive quick coupler	11
1.2	Technical Data: ProDrive Standard/Mini with light (ProDrive quick coupler)	7	2.2	Handpiece (Fixed connection)	12
1.3	Technical Data: ProDrive Standard/Mini without light (Fixed connection)	8	2.3	ProDrive Burs	13
1.4	Startup of ProDrive Standard/Mini	9	2.4	Setting the Cooling Spray	16
1.5	Notes on Safe Operation	9	3	Care and maintenance	17
1.6	Service life of ProDrive instruments	10	3.1	Care with Spray Lubricant (ProDrive quick coupler)	17

3.2 Care with Spray Lubricant (Fixed connection)	19
3.3 Care with Drop Oil Lubricant	21
3.4 Cleaning the Push-button Chuck	23
3.5 Cleaning & Disinfection	24
3.6 Sterilization	27
3.7 Cleaning The Cooling Spray Jets	29
4 Disposal	30

INDICATIONS FOR USE:

To be used intraorally by trained dental professionals for drilling and preparation of dental cavities for restoration, such as fillings.

CONTRAINDICATIONS:

Do not use this handpiece within the oral cavity when a soft tissue wound is present. Air pressure may force septic materials into these tissues.

WARNING: This product is intended for use by dental professionals only. Before operating handpiece, carefully read and follow these instructions and save them for future reference. Observe all cautions and warnings, including: Handpieces shall be heat-sterilized between patients to prevent cross-contamination which could result in serious illness or

even death from infectious organisms such as HIV (the virus that causes AIDS) or hepatitis B.

- Suitable hygienic measures must be taken to prevent cross contamination among patients, users and other persons.
- Do not cut with the handpiece unless you are sure the bur is secure. A loose or improperly extended bur could eject from the head or break and cause personal injury.

WARNING: ProDrive highspeed handpieces should not be rebuilt. Replacing bearings only could result in unexpected failure and patient injury due to different wear factors of new vs. old components.

1 Technical Description

1.1 ProDrive Standard and Mini

in compliance with ISO 7785 - 1 / ISO 9168

The ProDrive Standard/Mini highspeed handpieces and quick couplings comply with the relevant regulations, in accordance with the state of the art.


According to these regulations, such products may be used only for the described applications by a competent user and in compliance with the applicable hygiene standards, the relevant accident prevention measures and these operating instructions.

According to these regulations, the user is obliged to use only flawless working materials, solely for their intended purpose and to protect himself, the patient and any third parties against danger.

The ProDrive Standard/Mini highspeed handpieces and quick couplings are part of the technical equipment used in dental practices.

The quick coupling is used to supply media and couple the supply hose and is conveniently rotatable.


The ProDrive Standard/Mini highspeed handpieces are used to power dental instruments with a triangular drive shaft (\varnothing 1.6 mm).

ProDrive Standard**ProDrive Mini**

1.2 Technical Data: ProDrive Standard/Mini with light (ProDrive quick coupler)

2 different head sizes: Standard and Mini
3 different spray options: 1-jet, 2-jets, 3-jets,
delivered: 3-jets **with** illumination.

- Driving air pressure: 39 psi \pm 1.5 psi (2.7 bar \pm 0.1 bar) when flowing, dry and clean.
- Air consumption: 1.52 \pm 0.14 Cu ft/min (43 \pm 4 NI/min)
- Idling speed:
ProDrive Standard: 400,000 rpm \pm 10%
ProDrive Mini: 420,000 rpm \pm 10%.
- Prevention of suction into head housing.
- Anti-retraction valve.

ProDrive Standard**ProDrive Mini**

1.3 Technical Data: ProDrive Standard/Mini without light (Fixed connection)

2 different head sizes: Standard and Mini
 3 different spray options: 1-jet, 2-jets, 3-jets,
 delivered: 3-jets **without** illumination.

- Driving air pressure:
 Fixed connection: 32 psi \pm 1.5 psi (2.2 bar \pm 0.1 bar)
 when flowing, dry and clean.
- Air consumption: 1.52 \pm 0.14 Cu ft/min (43 \pm 4 NI/min)
- Idling speed:
 ProDrive Standard: 400,000 rpm \pm 10%
 ProDrive Mini: 420,000 rpm \pm 10%.
- Prevention of suction into head housing.
- Anti-retraction valve.

1.4 Startup of ProDrive Standard/Mini

Operating Pressure Parameters (flowing air)

- Driving air pressure: Handpiece with quick coupling:
39 psi \pm 1.5 psi (2.7 bar \pm 0.1 bar)
Driving air pressure: Handpiece with fixed connection:
32 psi \pm 1.5 psi (2.2 bar \pm 0.1 bar)
- Air exhaust: < 2 psi (< 0.15 bar)
- Spray air: Handpiece with quick coupling: 39 psi (2.7 bar)
Handpiece with fixed connection: 32 psi (2.2 bar)
- Spray water: 29 psi (2.0 bar)
- Recommended proportion of water in spray:
at least 50 ml/min.

The specified pressure values can be checked by a service engineer by inserting a testing device between the hose and the quick coupling.

1.5 Notes on Safe Operation

Sterilize the handpiece prior to initial operation.

Following longer **periods of disuse** the highspeed handpiece must be cleaned and cared for according to these instructions. During **prolonged preparations (>15 min.)** the highspeed handpiece should be treated only with spray lubricant.

Always store in a dry, heated room.

Discontinue use in case of irregular **running noises** or excessive vibrations and external signs of damage.

In case of noticeable malfunctions or damage, stop using the instrument immediately and contact your local dental depot or the manufacturer. If you have any questions, please also contact your dental depot or the manufacturer.

Use only approved or original replacement and accessory parts. Prior to the beginning of a workday and before each patient purge the water paths and air channels for 20 – 30s.

The logo for ProDrive, featuring the word "PRO" in blue, a stylized orange and blue circular icon, and the word "DRIVE" in blue.

1.6 Service life of ProDrive instruments


When used as intended:

- Non-moving parts of ProDrive instruments have a typical service life of 5 years*
- Moving parts of ProDrive instruments have a typical service life of 1 year*

Elastomers, e.g. O-rings, must be replaced depending on wear as described in the operating instructions.

For safety reasons, however, we recommend an annual check of the instrument holder system.

*No warranty claim can be inferred here, as wear may occur earlier or later than indicated above depending on use, frequency of sterilization and frequency of maintenance.

2 Operation

2.1 Handpiece with ProDrive quick coupler

ProDrive Standard and Mini

Snap on handpiece (1.) until it **audibly** engages.

Remove handpiece (2.) by holding quick coupling and hose nut securely while twisting the handpiece slightly and pulling it off toward the front.

⚠ CAUTION! *Do not remove handpiece while highspeed handpiece is running.
Do not hold by or pull on supply hose.*

2.2 Handpiece (Fixed connection)

ProDrive Standard and Mini


Connection to the supply hose


Attach the handpiece (1) to the supply hose, observing the various tube diameters.

Hold the handpiece (1) and tighten the cap nut (2) firmly by hand.

The sealing disk (3) between the supply hose and the handpiece can be replaced as required.

Note the **position of the bulging side (4)** when fitting the sealing disk (3).


2.3 ProDrive Burs

⚠ CAUTION! Use only flawless bits and diamond burs. These instruments should also occasionally be checked for sharpness and concentricity.

⚠ CAUTION: Make sure bur is secure before operating handpiece. Please refer to the manufacturer's specifications for the proper application and operational data of the various instruments.


- Total length up to max. 25 mm (ProDrive Standard)
- Total length up to max. 21 mm (ProDrive Mini)

ProDrive Burs

⚠ CAUTION! *You must use ProDrive burs in your ProDrive handpiece. Traditional, friction grip burs will not fit into the ProDrive Turbine. Attempting to forcibly insert friction grip burs may damage the Turbine and void your warranty.*

Insertion: When inserting the bur you will be required to depress the push-button on the handpiece as you rotate the bur so that it can be fully engaged. While inserting the bur, if you meet any resistance, simply rotate and pump the bur up and down as you insert it and it will fully engage.


Bur Lengthening

⚠ CAUTION! Make sure you do not accidentally actuate the push button by pulling back the cheek while the high-speed handpiece is running - danger of burning the patient's oral mucosa!

To lengthen the bur, first fully engage the bur as per the above instructions for bur insertion, then depress the push-button on the handpiece and slide the bur out until it is engaged into the upper bur safety groove. You will notice that the laser lengthening etch mark will be exposed when the bur is fully lengthened.

Removal:

⚠ CAUTION! After the highspeed handpiece comes to a stop press the push button and pull out the bur.

2.4 Setting the Cooling Spray

Set for the 6-Pin Coupler

using the control ring (4), while holding the hose nut (2) securely.


The **maximum water flow** is set when the control ring noticeably engages.

For the **5-Pin Coupler and the handpiece with fixed connection** the cooling spray must be set on the unit.

⚠ CAUTION! *A sufficient amount of coolant must always be available during preparation (proportion of water > 50ml/min).*

Deficient cooling will result in overheating of the preparation spot and damage to the tooth!

For cleaning the cooling spray jets see page 29.

3 Care and maintenance


3.1 Care with Spray Lubricant (ProDrive quick coupler)

Apply spray lubricant **prior to** every sterilization and **after** every thermal disinfection!

With prolonged use in one procedure, the handpiece should be lubricated at **15 minute intervals**.

- Remove the drilling instrument.
- Pull handpiece off of quick coupling.
- Push appropriate spray attachment (1) onto nozzle of spray can as shown.

i NOTE *The O-rings (2) must be undamaged.*


- Grasp the handpiece with a cloth and push it onto the spray attachment **up to the stop** and **hold it tight**.
- Spray for **1 - 2 seconds**.

If the fluid flowing out of the head is still dirty, spray again. Let the highspeed handpiece run briefly in between to expel excess spray lubricant.

- Wipe off emerging oil **when the highspeed handpiece is not running**.
- Before using the highspeed handpiece let it run briefly and wipe off any excess oil **when the handpiece is not running**.

- 1 Spray attachment
- 2 O-rings


3.2 Care with Spray Lubricant (Fixed connection)

Apply spray lubricant at least once every noon and evening, **prior to** every sterilization and **after** every thermal disinfection!

During **longer preparations (>15 min.)** the highspeed handpiece must be treated with spray lubricant in between.


- Remove the drilling instrument.
- Unscrew and pull off the handpiece of the supply hose.
- Make sure that the black nozzle adapter (1) is fitted in the spray can nozzle.

- **Hold** the handpiece using a cloth and insert the spray can nozzle in the driving air tube.
- Spray for **1 - 2 seconds**.


If the fluid flowing out of the head is still dirty, spray again. Let the turbine run briefly in between to expel excess spray lubricant.

- Wipe off emerging oil **when the handpiece is not running**.
- Before using the highspeed handpiece let it run briefly and wipe off any excess oil **when the handpiece is not running**.


3.3 Care with Drop Oil Lubricant

- Wipe handpiece to clean exterior of debris using alcohol.
- Insert appropriate coupler to handpiece.
- **Place 2 drops** of oil lubricant into the drive air tube.

- **Place 1 Drop of oil** into spindle opening.
- Insert ProDrive test bur to the fully seated position.
- Connect handpiece to the air station and operate at manufacturer's recommended operating specifications, for 30 seconds to purge all excess oil from the handpiece.
- Wipe off excess oil from handpiece with a wipe cloth. Remove ProDrive test bur from the handpiece.


3.4 Cleaning the Push-button Chuck

In order to eliminate deposits and preserve the function of the chucking system, the collet chuck should be cleaned only by **spray lubricant** at least once a week.

Attach highspeed handpiece head with chuck securely to nozzle of spray can and spray collet chuck for **1 - 2 seconds**.

Always hold the spray can upright!

Wipe off emerging oil with a lint-free cloth.

3.5 Cleaning & Disinfection

This procedure recommends the steps to follow between each patient use to prevent cross-contamination and ensure continued trouble-free operation of your ProDrive handpiece. Manual cleaning must always be combined with disinfection.

⚠ CAUTION! *The sterilization steps described in section 3.6 must be followed before every patient use whether or not a user has cleaned and disinfected their ProDrive handpiece.*



Cleaning the Handpiece

The highspeed handpiece can be brushed off under running water (of at least drinking water quality).

Cleaning Light Guide Surfaces

Blow surfaces (1) clean with compressed air to prevent scratching, dust or dirt particles etc.

Wipe off surfaces using a cotton pad or a soft cloth with alcohol.



Disinfection

The ProDrive Standard/Mini highspeed handpieces can be disinfected with steam at temperatures of up to 105° C (221°F), (3 min. holding time).

The ProDrive Standard/Mini highspeed handpieces are completely washable in the thermal disinfector at temperatures of up to 95°C.

Treat with spray lubricant afterwards. You must choose only instrument disinfectants whose bactericidal, fungicidal and virucidal properties have been tested and properly approved.

Spray and wipe disinfection of quick couplings and **ProDrive Standard/Mini highspeed handpiece** using chemical disinfectants is permissible only on the exterior surface!

You can use for instance:

CAVICIDE® or **CAVIWIPES™**

⚠ CAUTION! *Never immerse in disinfectant solutions!
Do not thermally disinfect quick couplings!*

277°F

3.6 Sterilization

The ProDrive Standard/Mini highspeed handpieces may be sterilized **only in autoclaves** with saturated water vapor at 136 °C (277 °F), at 220 kPa (2.2 bar, 31.9 psi) for 10 minutes.


Do not sterilize quick coupling!

⚠ CAUTION! *Do not exceed 284°F (140°C) during the drying phase!*


Approved for sterilization are steam sterilizers that at least use fractionated vacuum and are suitable for the sterilization of highspeed handpieces.

Before sterilizing:

Remove bur, Clean and disinfect highspeed handpiece (page 24), Apply spray lubricant to highspeed handpiece for protection.

**ProDrive
6-Pin Coupler**

**ProDrive
5-Pin Coupler**


277°F


⚠ CAUTION! Do not hasten cooling by submerging in water or placing the handpiece under running water. This will permanently damage your handpiece.

i NOTE Sterilize again once the storage period has elapsed.

**ProDrive
6-Pin Coupler**

**ProDrive
5-Pin Coupler**


3.7 Cleaning The Cooling Spray Jets

If calcareous tap water is used, the cooling spray jets may be partially or fully clogged due to calcification.

⚠ CAUTION! *A sufficient quantity of coolant must always be available during preparation (proportion of water >50 ml/min). Deficient cooling will cause overheating of the preparation spot, resulting in damage to the tooth.*

We therefore recommend that you **regularly** clean the jet openings **carefully** with the cleaning wire (1).

Then let the highspeed handpiece run briefly with cooling spray.

4 Disposal

To the best of our current knowledge, the product does not contain any substances which are hazardous to the environment.

The product must be sterilized before disposal.

Please observe the applicable national regulations regarding disposal.

ProDrive Standard/Mini

Turbines dentaires

PRODRIVE
SYSTEMS

Handpiece Performance Through Innovation

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ProDrive Systems Inc.
ProDrive Systems Ltd.